Test report

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1. Test process:**试验过程****1.1. Chromatographic conditions:****色谱条件**

Chromatographic column: 色谱柱	Ultimate ALK-C18 (4.6×250mm,5μm)
Mobile phase: 流动相	6.8 g/L potassium dihydrogen phosphate (pH adjusted to 6.0 with potassium hydroxide solution)-methanol-acetonitrile (65:15:20); 6.8g/L 磷酸二氢钾(氢氧化钾溶液调 pH 至 6.0)-甲醇-乙腈(65:15:20);

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1.2. Solution preparation (溶液配制)

1.2.1. Preparation of mobile phase (流动相配制):

Potassium hydroxide solution (5%): Weigh 5 g of potassium hydroxide, add 100 ml of purified water to dissolve it, and the solution is obtained.

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Mobile Phase : Measure out 650 ml of 6.8 g/L potassium dihydrogen phosphate solution (weigh 6.8 g of potassium dihydrogen phosphate, dissolve it in 1000 ml of pure water, ensure complete dissolution, adjust the pH to 6.0 with potassium hydroxide solution, and the solution is obtained), 150 ml of methanol, and 200 ml of acetonitrile. Mix them thoroughly, then filter and sonicate to obtain the solution.

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Diluent 2:Mobile Phase

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1.2.3. Test solution (测试溶液):

Stock solution: Weigh appropriate amounts of Compound D, Compound U, Compound T, Compound J, Compound 1, Compound N, and Compound K separately, dissolve them in methanol/methanol-water (60:40), and dilute to prepare a solution containing approximately 1 mg per mL.

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Test solution: Weigh an appropriate amount of Compound 3, separately pipette appropriate amounts of each of the above compound stock solutions, add mobile phase to dissolve and dilute to prepare a mixed solution containing approximately 1 mg of Compound 3, 10 µg each of Compound D, Compound U, Compound J, Compound K, and Compound N, and 20 µg each of Compound 1 and Compound T per mL.

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2. Chromatogram and data (谱图和数据)

(1) Analysis spectrum of Test solution (测试溶液分析图谱):


Figure 1: Acquired Results of the Test Solution at a Wavelength of 254 nm
图 1 测试溶液在 254nm 波长下采集结果


Figure 2: Acquired Results of the Test Solution at a Wavelength of 288 nm

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3. Conclusion (结论)

Using the Welch Materials Ultimate ALK-C18 chromatographic column (4.6×250 mm, 5 μm) under these chromatographic conditions, Compound 3 and its related substances can be effectively separated, and the separation between Compound 3 and D, U, and T all meet the requirements.

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